

AR@ptix
Switzerland

GENERAL PRODUCTS CATALOG



Table of Contents

ARCOPTIX S.A.	1
PRODUCTS LINE	2
SOFTWARE	3
A MATTER OF RESOLUTION	4
FREQUENTLY ASKED QUESTIONS	5
BENCHTOP INSTRUMENTS	8
ARCOPTIX FT-MIR ROCKET	9
ARCOPTIX FT-NIR ROCKET	11
ARCOPTIX FIBER COUPLED FT-IR.....	13
ARCOPTIX VIS-NIR-FIB SPECTROMETER.....	15
ARCOPTIX VIS-NIR-DR SPECTROMETER	17
OEM INSTRUMENTS	19
ARCOPTIX OEM010	20
ARCOPTIX OEM011	22
ACCESSORIES	24
ARCOPTIX COMMON PATH MODULE	25
ARCOPTIX LIQUID CELL MODULE	26
ARCOPTIX ARCLIGHT-MIR & ARCLIGHT-NIR.....	28
ARCOPTIX ARCSPIHERE-NIR	29
ARCOPTIX ARCSPIHERE-MIR	30
BATTERY OPERATED PORTABLE FTNIR SETUP	31
FIBERED REFLECTION PROBES.....	32
OPTICAL FIBERS.....	33
ARCOPTIX CUV	34
ARCOPTIX DRS98.....	35


ARCOPTIX S.A.

Who are we ?

ARCOptix S.A. is a swiss company active in the field of optical spectroscopy. Founded in 2006, our facilities are spread between our sales offices in Neuchâtel, located on the shore of lake Neuchâtel, and our production line in Tramelan, which lies in the beautiful region of the Jura mountains. Our core line of business is Fourier-transform infrared (FT-IR) spectroscopy, which can address a wide range of applications in the near-infrared (NIR) and mid-infrared (MIR) regions of the electromagnetic spectrum. With hundreds of instruments operating worldwide for more than a decade, our teams have an extensive expertise in the manufacturing and qualification of benchtop and original equipment manufacturer (OEM) FT-IR spectrometers.

What do we do ?

FT-IR spectrometers analyze the spectrum of light thanks to the scanning of an interferometer that lies at the core of the instrument. Our instruments rely on a permanently aligned interferometer, which consists in a pair of reflective corner-cubes mounted on a common swinging arm. This Y-shaped pendulum rotates on a wear-free flexure pivot, ensuring a long-lived operation by avoiding mechanical wear.



Dual corner cube mounted on a common swinging arm



Our OEM010 series featuring an internal light source, interferometer and detector.

This robust design paired with high quality, carefully selected optical components allow us to produce extremely compact, reliable and accurate analyzers with competitive performances with respect to other solutions found on the market of Fourier-transform spectroscopy. Besides standard products, we have a fast-prototyping capability that make possible the quick implementation and testing of custom configurations tailored to your specific application.

PRODUCTS LINE

**Diffuse reflection
 VIS-NIR**
 Page : 17



Fibered VIS-NIR
 Page : 15



VIS-NIR Units




NIR FT-IR

Desktop FT-IR
 Page : 11



NIR

MIR



Desktop FT-IR
 Page : 9



**Fibered FT-IR w/
 built-in light source**
 Page : 13




OEM FT-IR
 Pages : 20, 22

SOFTWARE

Our software – your application


At ARCOptix, we fully appreciate and value the multiple benefits that a dedicated, performant and reliable software can bring to your application. Automatic data collection, parameters changes, status diagnosis and many other essential tasks should be implemented as simply and as efficiently as possible in order to get the most out of your spectrometer. This philosophy led to the development of a multi-threading, cross-platform and versatile software application, the ARCOptix digital acquisition system or AoDAQ.



The AoDAQ simultaneously takes care of:

1. Handling communication with the FT-IR via USB
2. Processing raw signals to deliver a spectrum
3. Running a TCP Ethernet server

The AoDAQ can be installed on all sorts of computers, from desktop machines to embedded, low-power single board computers. Thanks to the hosting of a TCP server, the instrument data and parameters can be accessed locally and/or remotely. All communication with the instrument eventually reduces to a set of TCP/IP commands that allow to quickly acquire data, adjust parameters, monitor the instrument status etc. using the programming environment of your choice.



A MATTER OF RESOLUTION

Which resolution for your application ?

The resolution of an optical spectrometer is without doubt one of its most crucial parameters. Higher resolving power means better discrimination of closely positioned spectral features, which intuitively speaking should enable retrieving more detailed information from a given measurement. While this is true to some extent, some applications do not benefit further from a better resolution beyond a certain point. Besides, the resolution of an FT-IR cannot be improved without trading other parameters such as acquisition time or signal-to-noise ratio, thus degrading other important signal metrics. Overall, what matters most is to operate at the *right* resolution, not at the sharpest available one.

At ARCOptix, our instruments are classified in two distinct categories. Standard resolution instruments offer tunable resolutions up to 2cm^{-1} , which is ideal for the study of condensed matter in general, as it matches the width of the spectral features observed in solids and liquids in the mid and near infrared. High resolution instruments can reach resolution up to 0.5cm^{-1} , which enable resolving the finer structure of individual molecules in a gas mixture or can provide more insight on the emission of a given source such as a laser.

Standard resolution (up to 2cm^{-1}),
recommended for:




Solid measurements



Liquid measurements

High resolution (up to 0.5cm^{-1}),
recommended for:



Gas spectroscopy




Laser characterization

FREQUENTLY ASKED QUESTIONS

What is the equivalent wavelength resolution of my FT-IR ?

Owing to its operating principle, a FT-IR delivers a spectrum uniformly sampled over a given spectral range in wavenumbers (ν), with units of cm^{-1} . The wavenumber is simply defined as the inverse of the wavelength (λ). The resolution of an FT-IR is constant in wavenumber ($\Delta\nu$), but varies with wavelength ($\Delta\lambda$) due to the inverse relationship between these two units. The equivalence is given by the following rule:

$$\Delta\lambda = \frac{\Delta\nu}{\nu^2}$$




Resolution in [nm] as a function of wavelength in FT-IR systems

What is the typical emission of my ARCLIGHT-MIR/NIR ?


The ARCOptix ARCLIGHT lamps are thermal sources, which behave to a good approximation like blackbodies at a temperature of 2'850 °K (ARCLIGHT-NIR) and 1'450 °K (ARCLIGHT-MIR), the spectral radiances of which are illustrated below. The typical total power emitted by both version of the ARCOptix ARCLIGHT source are shown in the following table.

	ARCLIGHT-NIR	ARCLIGHT-MIR
w/o fiber coupler	360 mW	180 mW
w/ fiber coupler	110 mW	40 mW
w/ optical fiber	20 mW	4 mW

Total optical power of the ARCLIGHT in different configurations




Typical radiance of the ARLIGHT-NIR



Typical radiance of the ARCLIGHT-MIR

How to choose the best detector for my application ?


The choice of a detector is a most sensitive one as it contributes to a large extent to the overall quality of your measurement. ARCOptix offers several types and models of detectors, namely extended InGaAs detectors (2 stage thermo-electric cooling), HgCdTe (MCT) detectors with either thermoelectric cooling (4 stages) or liquid nitrogen cooling as well as pyroelectric (DLATGS) detectors. Selection of a given detector is application driven, and depends on the spectral signatures that you are willing to analyze. In order to deliver a faire comparison between various types of photodetectors (or thermal detectors), the most widespread metric is the specific detectivity (labelled D^*), which is given in the figure below for all of our available detectors.




Specific detectivity (D^*) of ARCOptix detectors

What configuration should I select for my application?

ARCOptix offers two distinct classes of spectrometers: FTIRs without a built-in light source, such as the **FTMIR-L1-120-4TE**, and FTIRs which integrate a built-in emitter, such as the **FTMIR-FC-120-4TE**. Source-free instruments can be used to characterize both samples (relative measurements) and optical sources (absolute measurements). FTIR with built-in sources can only be used to perform relative measurements, as light is modulated by the interferometer before being output. This configuration has the advantage of being robust to signal modification induced by stray-light or background thermal radiation-



Light routing in the FTMIR-L1 -120-4TE (left) and the FTMIR-FC-120-4TE (right)




BENCHTOP INSTRUMENTS

Are you a scientist, engineer, technician or scholar looking for a compact, reliable and efficient FT-IR solution ? Our portfolio of benchtop FT-IR instruments can target a wide range of applications in the near-infrared or mid-infrared, such as liquid analysis, solid sample characterization, gas spectroscopy, diffuse reflectance measurements and many more. Our benchtop instruments are designed for ease of use and flexibility, and come with a large choice of sources and accessories to help you get the information that you need.

- ❖ [ARCOPTIX FT-MIR ROCKET: Page 9](#)
- ❖ [ARCOPTIX FT-NIR ROCKET: Page 11](#)
- ❖ [ARCOPTIX FIBER COUPLED FT-IR: Page 13](#)
- ❖ [ARCOPTIX VIS-NIR-FIB SPECTROMETER: Page 15](#)
- ❖ [ARCOPTIX VIS-NIR-FIB SPECTROMETER: Page 17](#)

ARCOPTIX FT-MIR ROCKET



The ARCOptix FT-MIR Rocket is a highly performant, compact and reliable spectrometer that is ideal for various applications in the mid-infrared. The concentration levels of CO₂ or H₂O in the interferometer volume are conveniently minimized thanks to a homemade, replaceable dessicant capsule. Thanks to its permanently aligned interferometer and solid-state reference laser, the FT-MIR Rocket offers excellent stability in both intensity and wavelength scales.

With four available spectral ranges and adjustable spectral resolution down to 2cm⁻¹ (0.5cm⁻¹ on request), the FT-MIR Rocket is a highly flexible instrument that can be tailored to your application. Designed for convenience and ease-of-use, our FT-MIR spectrometer is readily operational with our ARCSpectroRocket software using a standard USB 2.0 connection.

Applications

- *Mid-IR Optical Spectrum Analyzer (OSA) for MIR Lasers & LEDs*
- *Liquid, thin-film or gas measurement*
- *Material identification and quantification in various fields such as geology, food and beverage industry, ...*

Features

- **4 spectral ranges :**
 - 2-6 μm (TEC-MCT)
 - 1.5-8.5 μm (TEC-MCT)
 - 2-12 μm (TEC-MCT)
 - 2-16 μm (LN2C-MCT or DLATGS)
- **Dynamically adjustable resolution:**
 - 8cm⁻¹
 - 4cm⁻¹
 - 2cm⁻¹
 - 0.5cm⁻¹ (on request)
- **Compact design:**
18cm X 16cm X 8 cm
- **Wear free moving parts for extended lifetime**
- **No purging of the interferometer required**
- **Removable fiber adapter**
- **Temperature controlled reference laser**
- **Low power consumption**
- **USB 2.0 connection**

Specifications

Product code	FTMIR-L1-060-4TE	FTMIR-L1-085-4TE	FTMIR-L1-120-4TE	FTMIR-L1-160-LN2	FTMIR-L1-160-DLA
Beamsplitter material	CaF ₂		ZnSe		
Spectral Range [cm ⁻¹]	5'000 – 1'660	6600-1'200	5'000 - 830	5'000 - 650	
Spectral Range [μm]	2-6	1.5-8.5	2-12	2-16	
Detector Type	MCT (4-TE cooled)			MCT (LN2 cooled)	DLATGS
Detector Peak D* [cm Hz ^{1/2} W ⁻¹]	>1x10 ¹¹	>8x10 ⁹	>4x10 ⁹	>5x10 ¹⁰	>2.5x10 ⁸
Signal-to-noise ratio	> 80'000:1 ⁱ	> 40'000:1 ⁱ	> 40'000:1 ⁱ	>70'000:1 ⁱ	>8000:1 ⁱ
Removable fiber-optic coupler	Lensed (CaF ₂ fiber coupler)		Reflective fiber coupler (90° off-axis parabolic mirror)		
Recommended fiber	IFG (1-6μm)	IFG (1-6μm) or PIR (3-18μm)	PIR (polycrystalline) fibers, 3-18μm		
Fibered interface	Fiber core up to Ø 900μm, NA=0.3, SMA 905 connector				
Free-space interface	Ø 12.7mm collimated (max ~30mrad half angle)				
Interferometer type	Permanently aligned, double retro-reflector design				
Resolution (unapodized) [cm ⁻¹]	0.5 ⁱⁱ , 2, 4, 8 (user selectable)				
Wavenumber repeatability	<10 PPM				
Scan frequency	>4 Hz @ 4cm ⁻¹				>0.4Hz @ 4cm ⁻¹
Internal reference laser	Temperature-stabilized solid-state laser @850nm				
A/D Converter	24 bit				
Operating temperature	10°C-40°C				
Power requirement	12V / 10W max			12V / 6W max	
Communication Interface	USB 2.0				
Software Interface	Windows 10/11 API for controlling the instrument via our DLL				
Dimensions	180mm x 160mm x 80mm (without Dewar)				
Weight	1800 g (without Dewar)				

ⁱ Measured with a silicon carbide (SiC) source (~1550°K) with f=18mm reflector directly shining into the free-space input port, 60s measurement, around peak sensitivity wavelength, 4cm⁻¹, Norton-Beer weak apodization.

ⁱⁱ Available on request only, please contact us at info@arcoptix.com for details

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE. Please contact info@arcoptix.com for more information.

ARCOPTIX FT-NIR ROCKET



The ARCOptix FT-NIR Rocket is a highly performant, compact and reliable spectrometer that is ideal for various applications in the near-infrared spectrum. Its fibered port makes it directly compatible with fibered accessories such as reflection probes or cuvette holders to analyze e.g. liquids in transmission.

Thanks to its permanently aligned interferometer and solid-state reference laser, the FT-NIR Rocket offers excellent stability in both intensity and wavelength scales. With four available spectral ranges and adjustable spectral resolution down to 2cm^{-1} , the FT-NIR Rocket is a highly flexible instrument that can be tailored to your application. Designed for convenience and ease-of-use, our FT-NIR spectrometer is readily operational with our ARCSpectroRocket software using a standard USB 2.0 connection.

Applications

- *Transmission, diffuse reflectance*
- *Light source measurement (NIR Lasers, LED, Solar,...)*
- *Material identification and quantification in various fields such as geology, food and beverage industry, drug & medical diagnostics, etc.*

Features

- **High D* 2-TE cooled InGaAs photodetector**
- **Extended range available with dual detector**
- **Dynamically adjustable resolution:**
 - 8cm^{-1}
 - 4cm^{-1}
 - 2cm^{-1}
- **Compact design:**
18cm X 16cm X 8 cm
- **Wear free moving parts for extended lifetime**
- **No purging of the interferometer required**
- **Low power consumption**
- **Compatible with various fibered accessories**
- **Temperature controlled reference laser**
- **USB 2.0 connection**

Specifications

Product code	FTNIR-L1-025-2TE	FTNIR-L1-060-EXT
Spectral Range [cm⁻¹]	11'000-4'000	11'000-1'660
Spectral Range [μm]	0.9-2.5	0.9-6.0
Detector type	Extended type InGaAs 2TEC	Extended type InGaAs 2TEC & MCT 4TEC
Detector Peak D* [cm Hz^{1/2}W⁻¹]	>2x10 ¹¹	>2x10 ¹¹ (0.9-2.5μm) >1x10 ¹¹ (2.5-6.0μm)
Signal-to-noise ratio (SNR)	>100'000:1 ⁱ	>100'000:1 ⁱ (0.9-2.5μm) >80'000:1 ⁱ (2.5-6.0μm)
Removable fiber-optic coupler	Lensed (CaF2 fiber coupler)	
Recommended fiber	Low-OH silica Multi-mode fiber with 200μm or 600μm core	-
Fibered interface	Fiber core up to Ø 0.6mm, NA=0.25, SMA 905 connector	
Free-space interface	Ø 12.7mm collimated (max ~30mrad half angle)	
Interferometer type	Permanently aligned with dual retro-reflector	
Resolution (unapodized) [cm⁻¹]	2, 4, 8 (user selectable)	
Wave-number repeatability	<20 PPM	
Scan frequency	>4 Hz @ 4cm ⁻¹	
Internal reference laser	Temperature controlled solid-state @795nm	
A/D Converter	24 bit	
Amplifier	4 gain levels low noise trans-impedance amplifier	
Operating temperature/ humidity	5°C - 40°C / non condensing	
Power requirements	12V / 10W max	12V / 20W max
Communication interface	USB 2.0	
Software Interface	Windows 10/11 API for controlling the instrument via our DLL	
Dimensions	180mm x 160mm x 80mm	
Weight	1800 g	

ⁱMeasured with a 20W halogen lamp in transmission mode, 60s measurement, around peak sensitivity wavelength, Norton-Beer weak apodization, linearly corrected baseline, resolution setting 4 cm⁻¹

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE. Please contact info@arcoptix.com for more information.

ARCOPTIX FIBER COUPLED FT-IR



The ARCOptix FT-IR spectrometer is an all-fibered alternative to our “Rocket” configurations. It features the same components as the FT-MIR Rocket, with the addition of an internally mounted optical source that enables to modulate light before coupling it to the output port of the device. This architecture is thus more robust to external perturbations such as background thermal emission, as such parasitic light would not go through the interferometer before reaching the detector. Our FTIR-FC is ideal for fibered applications with the exception of diffuse reflectance probes, which benefit from the higher throughput provided by direct feeding from a standalone source such as our ARCLIGHT-MIR.



FTMIR-FC unit mounted with an MCT detector cooled by LN2

Features

- **Internal light source**
- **SMA905 fiber connectors**
- **Available with LN2 cooled MCT photodetector**
- **Robust to ambient light perturbation**
- **Dynamically adjustable resolution:**
 - 8cm^{-1}
 - 4cm^{-1}
 - 2cm^{-1}
- **Wear free moving parts for extended lifetime**
- **No purging of the interferometer required**
- **Temperature controlled reference laser**
- **Low power consumption**
- **USB 2.0 connection**

Specifications

Product code	FTMIR-FC-060-4TE	FTMIR-FC-120-4TE	FTMIR-FC-160-LN2
Beam-splitter material	CaF ₂	ZnSe	
Spectral Range [cm ⁻¹]	5'000-1'660	5'000 – 830	5'000 – 650
Spectral Range [µm]	2-6	2-12	2-16
Detector Type	MCT (4-TE cooled)		MCT (LN2 cooled)
Detector peak D* [cm Hz ^{1/2} W ⁻¹]	>1 x10 ¹¹	>4x10 ⁹	>5x10 ¹⁰
Signal-to-noise ratio	> 80'000:1 ⁱⁱ	> 40'000:1 ⁱⁱ	> 70'000:1 ⁱⁱ
Recommended fiber	IFG (1-6 µm)	PIR (3-18µm)	
Fibered interface	Fiber core up to Ø 900µm, NA=0.3, SMA 905 connector		
Internal reference laser	850nm		
Power requirement	40W @12VDC		30W @12VDC
Integrated light source	20W SiC globalar		
Interferometer type	Permanently aligned, double retro-reflector design		
Resolution (unapodized) [cm ⁻¹]	2, 4, 8 (user selectable)		
Wavenumber repeatability	<10PPM		
Scan frequency	>4 Hz @ 4cm ⁻¹		
A/D Converter	24 bit		
Operating temperature	10°C-40°C		
Communication Interface	USB 2.0		
Software Interface	Windows 10/11 API for controlling the instrument via our DLL		
Dimensions	180mm x 160mm x 80mm		
Weight	2200 g (excluding LN2 dewar)		

ⁱMeasured with a 20W halogen lamp in transmission mode, 60s measurement, around peak sensitivity wavelength, Norton-Beer weak apodization, linearly corrected baseline, resolution setting 4 cm⁻¹

ⁱⁱMeasured with a silicon carbide (SiC) source (~1550K) in transmission mode, 60s measurement, around peak sensitivity wavelength, 4cm⁻¹, Norton-Beer weak apodization.

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE. Please contact info@arcoptix.com for more information.

ARCOPTIX VIS-NIR-FIB SPECTROMETER



The ARCOptix VIS-NIR spectrometer is a broadband spectroscopy solution covering both the visible (down to 350 nm) as well as the near-infrared (up to 2500 nm). This all-fibered unit integrates our ARCOptix FT-NIR spectrometer as well as a multichannel grating spectrometer to probe light in the visible range. The two individual spectra are then processed and merged automatically in our software to deliver a single measurement of the VIS-NIR spectrum.

Featuring a high resolution (better than 1.5 nm) over its entire measuring range, the ARCOptix VIS-NIR spectrometer is a highly convenient, compact and user friendly instrument for various applications. The unit includes a dedicated Y-branch fiber bundle optimized for each spectral range in order to maximize the signal measured by each built-in instrument.

Applications

- *Transmission measurement*
- *Diffuse reflectance measurement*
- *Light source characterization*
- *Material identification*

Features

- **Grating spectrometer & FT-NIR combined**
- **Y-branch fiber bundle**
- **Lightweight and compact (22x18x8 cm)**
- **High resolution < 1.5nm over the 350-2500nm range**
- **Fast acquisition time (2 s)**
- **USB 2.0 connection**
- **Access to instrument individual data & combined information**
- **Compatible with various fibered accessories**

Specifications

Product code	VIS-NIR-FIB																				
Spectral range	350nm-2'500nm																				
Resolution	< 1.5 nm																				
Fiber inputs	2X SMA-905																				
Single acquisition scanning time [s]	2																				
Software interface	Windows 10/11																				
Power supply	12 VDC (power supply included)																				
Communication interface	USB 2.0																				
Dimensions [mm]	220x180x80																				
Weight [kg]	2.5																				
Operating temperature [°C]	5-35																				
SNR	>1'000:1 @600nm (VIS) >10'000:1 @2000nm (NIR) NOTICE: Due to poor sensitivity of Si and InGaAs detectors, the SNR is up to 10x lower in the overlap region from 900 to 1100 nm.																				
<table border="1"> <caption>Grating Efficiency Data (Estimated from Graph)</caption> <thead> <tr> <th>Wavelength [nm]</th> <th>Grating efficiency [%]</th> </tr> </thead> <tbody> <tr><td>350</td><td>58</td></tr> <tr><td>400</td><td>62</td></tr> <tr><td>450</td><td>68</td></tr> <tr><td>500</td><td>70</td></tr> <tr><td>600</td><td>62</td></tr> <tr><td>700</td><td>52</td></tr> <tr><td>800</td><td>42</td></tr> <tr><td>900</td><td>35</td></tr> <tr><td>1000</td><td>30</td></tr> </tbody> </table>		Wavelength [nm]	Grating efficiency [%]	350	58	400	62	450	68	500	70	600	62	700	52	800	42	900	35	1000	30
Wavelength [nm]	Grating efficiency [%]																				
350	58																				
400	62																				
450	68																				
500	70																				
600	62																				
700	52																				
800	42																				
900	35																				
1000	30																				
VIS Spectrometer																					
Technology	Array detector grating spectrometer																				
Spectral range	350nm-1'000nm																				
Resolution	~1.5 nm (25 µm slit)																				
NIR Spectrometer																					
Technology	Fourier-transform spectrometer, model FTNIR-L1-025-2TE																				
Spectral range	900nm-2'500nm (extended InGaAs detector, 2TE cooled)																				
Resolution	<1.5nm (2cm ⁻¹)																				

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE. Please contact info@arcoptix.com for more information.

ARCOPTIX VIS-NIR-DR SPECTROMETER



The ARCOptix VIS-NIR-DR spectrometer is a broadband instrument designed for diffuse reflectance measurements in the visible and near-infrared. The entrance port of the instrument leads to an integrating sphere with internal illumination, which ensures optimized illumination of your sample. The VIS-NIR-DR integrates a multichannel grating spectrometer for measuring in the visible as well as an FT-IR that covers the NIR region of the electromagnetic spectrum.

Designed for ease of use, our dedicated software merges the information delivered by the two instruments and outputs a single spectrum covering the 350-2500nm interval. This software also features a direct calculation of the total solar reflectance (TSR) value for reference measurements.

Applications

- *Diffuse reflectance measurement of paints, textiles, ink, cosmetics, plastics etc.*
- *Measurement of total solar reflectance (TSR) in accordance to ASTM norms G173 and E903*

Features

- **Grating spectrometer & FT-NIR combined**
- **Sample holder**
- **Internally illuminated integrating sphere**
- **High resolution < 5nm over the 350-2500nm range**
- **Fast acquisition time (5 s)**
- **USB 2.0 connection**
- **Direct calculation of TSR value**

Specifications

Product code	VIS-NIR-DR
Spectral range [nm]	350-2500
Resolution [nm]	< 5
Measurement geometry	Diffuse illumination, 8° viewing angle
Integrating sphere diameter [mm]	50
Integrating sphere port diameter [mm]	10
Illumination source	Halogen, 5W
Single acquisition scanning time [s]	< 5
SNR (single acquisition)	> 1000:1
Inter-instrument repeatability	±1% at 50% R
Repeatability on white reference	<0.2% variation at 50% R
Software interface	Windows 10/11
Operating power	20 W @ 12 VDC
Communication interface	USB 2.0
Dimensions [mm]	380x250x320
Weight [kg]	8
Operating temperature [°C]	10-40

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE. Please contact info@arcoptix.com for more information.



OEM INSTRUMENTS

Are you looking to develop an industrial application based on Fourier-Transform spectroscopy? Compact, robust and reliable, our product line of OEM modules has been specifically designed for easy integration in a complete measuring system. Besides our standard portfolio, our engineering department has the experience and manufacturing capability to offer custom, dedicated solutions that can be tailored to match your current needs and requirements.


❖ [ARCOPTIX OEM010: Page 20](#)

❖ [ARCOPTIX OEM011: Page 22](#)

ARCOPTIX OEM010



The ARCOptix OEM010 is our most compact FT-IR complete system. It features a built-in SiC globar heated to 1550 °K and regulated in temperature to ensure stable illumination. The light that is being output by the instrument is first modulated by our permanently aligned interferometer, making it robust against stray light that might eventually reach the detector. The OEM010 rugged casing also holds a highly sensitive TE-MCT detector which maximizes the response of the system.



Optical routing in the ARCOPTIX OEM010

Features

- Internal light source
- Compact, rugged unit
- Custom I/O coupling modules on request
- Robust to ambient light perturbation
- Dynamically adjustable resolution:
 - 8cm⁻¹
 - 4cm⁻¹
 - 2cm⁻¹
 - 0.5 cm⁻¹ (on request)
- Wear free moving parts for extended lifetime
- No purging of the interferometer required
- Temperature controlled reference laser
- Ethernet interface with embedded processing available

Specifications

Product code	FTIR-OEM010-060-4TE	FTIR-OEM010-085-4TE	FTIR-OEM010-120-4TE	FTIR-OEM010-160-LN2
Beam-splitter material	CaF ₂		ZnSe	
Spectral Range [cm ⁻¹]	5'000-1'660	6'600-1'200	5'000-830	5'000-650
Spectral Range [µm]	2-6	1.5-8.5	2-12	2-16
Detector Type	MCT (4-TEC)			MCT (LN2)
Detector Peak D* [cm Hz ^{1/2} W ⁻¹]	>1x10 ¹¹	>8x10 ⁹	>4x10 ⁹	>5x10 ¹⁰
Signal-to-noise ratio	> 80'000:1 ⁱ	>40'000:1 ⁱ		>70'000:1 ⁱ
Output beam	Ø 12.7mm collimated (max ~30mrad half angle)			
Interferometer type	Permanently aligned, double retro-reflector design			
Resolution (unapodized) [cm ⁻¹]	0.5 ⁱⁱ , 2, 4, 8 (user selectable)			
Wavenumber repeatability	<10 PPM			
Scan frequency	>4 Hz @ 4cm ⁻¹			
Reference laser wavelength [nm]	850			
A/D Converter	24 bit			
Operating T [°C]	10 to 40			
Storage T [°C]	-20 to 60			
Built-in light source	SiC globar (1550°K)			
Power requirement	12V / 40W			
Communication Interface	USB 2.0 / Ethernet with embedded processing			
Software Interface	Windows 10/11, Linux API for controlling the instrument via our DLL Single board computer available on request			
Dimensions [mm]	165x145x82 (without dewar)			
Weight [g]	2100 (without dewar)			


ⁱInput to output direct light routing, 60s measurement, around peak sensitivity wavelength, 4cm-1, Norton-Bier weak apodization.

ⁱⁱ Available on request only, please contact us at info@arcoptix.com for details
 SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE.

ARCOPTIX OEM011



The ARCOptix OEM011 is a flexible alternative to our OEM010 series. The main module features a built-in light source (SiC globar) regulated in temperature as well as our permanently aligned interferometer system. The TE-MCT detector has been moved to an external module which is ideal for configurations requiring a sampling system (short path gas cell, purged volume, etc.). Both modules are easily fixed on optical breadboards and can accommodate 30 mm cage system rods for rapid prototyping.



Optical routing in the ARCOPTIX OEM011

Features

- Internal light source
- Compact, rugged unit
- External detector module
- Robust to ambient light perturbation
- Dynamically adjustable resolution:
 - 16cm^{-1}
 - 8cm^{-1}
 - 4cm^{-1}
 - 2cm^{-1}
 - 0.5cm^{-1} (on request)
- Wear free moving parts for extended lifetime
- Compatible with 30 mm cage system rods
- Temperature controlled reference laser

Specifications

Product code	FTIR-OEM011-060-4TE	FTIR-OEM011-085-4TE	FTIR-OEM011-120-4TE	FTIR-OEM011-160-LN2
Beam-splitter material	CaF ₂		ZnSe	
Spectral Range [cm ⁻¹]	5'000-1'660	6'600-1'200	5'000-830	5'000-650
Spectral Range [µm]	2-6	1.5-8.5	2-12	2-16
Detector Type	MCT (4-TEC)			MCT (LN2)
Detector Peak D* [cm Hz ^{1/2} W ⁻¹]	>1x10 ¹¹	>8x10 ⁹	>4x10 ⁹	>5x10 ¹⁰
Signal-to-noise ratio	> 80'000:1 ⁱ	>40'000:1 ⁱ		>70'000:1 ⁱ
Output beam	Ø 12.7mm collimated (max ~30mrad half angle)			
Interferometer type	Permanently aligned, double retro-reflector design			
Resolution (unapodized) [cm ⁻¹]	0.5 ⁱⁱ , 2, 4, 8 (user selectable)			
Wavenumber repeatability	<10 PPM			
Scan frequency	>4 Hz @ 4cm ⁻¹			
Reference laser wavelength [nm]	850			
A/D Converter	24 bit			
Operating T [°C]	10 to 40			
Storage T [°C]	-20 to 60			
Built-in light source	SiC globar (1550°K)			
Power requirement	12V / 30W (interferometer), 12V / 10W (detector module)			
Communication Interface	USB 2.0 / Ethernet with embedded processing			
Software Interface	Windows 10/11, Linux API for controlling the instrument via our DLL Single board computer available on request			
Dimensions [mm]	165x145x82 (interferometer), 93x75x66 (detector module, without dewar)			
Weight [g]	2100 (interferometer), 400 (detector module, without dewar)			

ⁱ Input to output direct light routing, 60s measurement, around peak sensitivity wavelength, 4cm⁻¹, Norton-Bier weak apodization.

ⁱⁱ Available on request only, please contact us at info@arcoptix.com for details
 SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE.




ACCESSORIES

Accessories for FT-IR systems are a must have to deliver stable, accurate and repeatable measurements. Naturally matched to our spectrometer product line, our various accessories will facilitate the setup of your specific experimental configuration.

- ❖ **ARCOPTIX COMMON PATH MODULE: Page 25**
- ❖ **ARCOPTIX LIQUID CELL MODULE: Page 26**
- ❖ **ARCOPTIX ARCLIGHT-MIR & ARCLIGHT-NIR: Page 28**
- ❖ **ARCOPTIX ARCSPHERE-NIR: Page 29**
- ❖ **ARCOPTIX ARCSPHERE-MIR: Page 30**
- ❖ **BATTERY OPERATED PORTABLE FTNIR SETUP: Page 31**
- ❖ **OPTICAL FIBERS: Page 33**
- ❖ **ARCOPTIX CUV: Page 34**
- ❖ **ARCOPTIX DRS98: Page 35**

ARCOPTIX COMMON PATH MODULE



Common path module operating principle



Common path module mounted on the OEM010

The ARCOptix common path (OEM-PART-CMP) module is an input/output accessory that can be mounted on our OEM010 series. The common path module features a beamsplitter that allows to guide the modulated beam output by the instrument and the light returning from the sample along the same optical path. This accessory coupled to our parabolic reflector for diffuse measurement (OEM-PART-DRP) offers a high-throughput, free-space solution for measuring diffuse and specular reflection in the MID-IR over a broad spectral range, without being limited to the transparency regions of CIR and PIR optical fiber.



ARCOptix FTIR OEM010 module mounted with the OEM-PART-CMP and the OEM-PART-DRP for diffuse reflection measurement.





Illustration of the OEM-PART-DRP focusing power.

Specifications

Product code	OEM-PART-CMP	OEM-PART-DRP
Internal optics	ZnSe beamsplitter	Gold parabolic reflector
Output	Collimated	Focused
Working distance	up to 1m	10mm (to focus from edge)

ARCOPTIX LIQUID CELL MODULE



Assembly consisting of the FTIR-OEM010-120-4TE, OEM-PART-CMP and OEM-PART-LCX



Close up view of the OEM-PART-LCX



The ARCOptix liquid cell module (OEM-PART-LCX) has been originally developed for application in high performance liquid chromatography (HPLC). However, the application range is much broader - it can be used as a standalone transmission cell for liquid or even high pressure gases. The OEM-PART-LCX is a high pressure, chemically and heat resistant cell for liquids spectroscopy. The cell has been tested to function up to 210 Bars and 180°C. The cell's body is made of L316 stainless steel (Hastelloy available on request). The optical path length can be set in the range 0.5mm-3.0mm. The wetted surfaces are: gold, sapphire (as a standard), stainless steel (or Hastelloy) and carbon. The choice of these materials makes it particularly suitable for spectroscopy in the transmission windows of organic solvents. The cell can be installed on a tubular extender, and thus it can be inserted inside a thermostat (for example directly next to an HPLC cell) avoiding the re-mixing of the separated mixtures. The spectral range is determined by the used detector of the FTIR and by the window of the cell. The standard model is operating in the wavelength range of 2 to 6 microns (limited by the sapphire window). For other applications it can be extended to 2-12 microns.

Specifications

Product code	OEM-PART-LCX
Path length	0.5 to 3 mm (chosen and fixed)
Internal volume	< 15 µL
Total transmission	<90% (with air inside)
Internal temperature [°C]	20-180
Mirrors	gold coating
Window material	Sapphire
Transmission range [µm]	2-6 (Sapphire)
Liquid inlet/outlet	Waters compatible fitting with 1/16" capillaries
Power requirement	35 W @ 110-230 VAC or 12 VDC
Dimensions [mm]	73x25Ø (without extender)

Measurement example

The pictures below shows a dilution sequence of mineral oil (10'000 ppm to 78 ppm) in a 2,2,3,3 Tetrachlorohexafluorobutane solvent with an optical path of 1.5 mm.



ARCOPTIX ARCLIGHT-MIR & ARCLIGHT-NIR



ARCLIGHT-MIR with reflector fiber coupler



ARCLIGHT-NIR with lens fiber coupler

The Arcoptix ARCLIGHT-MIR & ARCLIGHT-NIR are versatile lamps optimized for the VIS-NIR (400 nm – 4,000 nm) or MIR (1-25 μ m) spectral ranges. The lamps feature an electronically stabilized power supply, a variable mechanical attenuator, and are supplied with a removable fiber coupler enabling free-space or fiber-coupled operation.

Specifications

Product code	ARCLIGHT-NIR	ARCLIGHT-MIR
Light emitting element	QTH (Halogen)	SiC global
Source temperature [K]	2'850	1'550
Spectral range [μ m]	0.4-4	1-25
Spectral range [cm^{-1}]	25'000-2'500	10'000-400
Rated lifetime [hours]	4'000	10'000
Removable fiber-optic coupler	NA = 0.25, CaF2 lens	NA = 0.3, gold-coated off-axis parabolic mirror
Power consumption [W]	20 (electrical)	
Power requirement	12V via OD 2.5mm/ID2.1mm power jack (100-240 VAC to 12 VDC converter included)	
Output beam diameter [mm]	12.7	
Operational temperature [°C]	5-40	
Attenuator	Manual Iris Diaphragm	
Dimensions [mm]	92x92x43 (w/o fiber coupler), 140x92x43 (w/ fiber coupler)	
Weight [kg]	0.4	

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
 Please contact info@arcoptix.com for more information.

ARCOPTIX ARCSPHERE-NIR



The ARCOptix ARCSPHERE-NIR is an integrating sphere featuring an internal halogen light source for diffuse reflection measurements in the VIS-NIR. The highly reflective sphere has a diameter of 50mm, a sample port of 10mm with a sapphire window, and includes an SMA fiber connector allowing a direct connection to the spectrometer. The internally illuminated design provides a much stronger signal (at least 50 times more) than spheres illuminated via fiber optic coupled light sources. The internal volume of the sphere is cleaned by a desiccant capsule allowing to reduce significantly the contribution from water vapor on the measurement.

Specifications

Product code	ARCSPHERE-NIR
Sphere internal diameter [mm]	50
Sample port diameter [mm]	10
Sample port window material	Sapphire
Optical fiber connector	SMA-905
Light emitting element	Low power halogen bulb (5W)
Source temperature [K]	~2'700
Rated life-time [hours]	4'000
Power consumption (7-12 VDC) [W]	10
External dimensions [mm]	70x70x90
Weight [g]	900

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
Please contact info@arcoptix.com for more information.

ARCOPTIX ARCSPHERE-MIR



The ARCOptix ARCSPHERE-MIR is an integrating sphere featuring an internal SiC globar for diffuse reflection measurements in the MIR. Light emitted by the source is prevented from illuminating the sample directly by the presence of a deflector (baffle). The highly reflective gold coated sphere has a diameter of 50mm and a sample port of 10mm with a BaF₂ window. The output of the sphere is directly injected in one of our FTIR, such as the **FTMIR-L1-120-4TE**, maximizing throughput. The internal volume of the sphere is cleaned by a desiccant capsule allowing to reduce significantly the contribution from water vapor and CO₂ on the measurement.

Specifications

Product code	ARCSPHERE-MIR
Sphere internal diameter [mm]	50
Sample port diameter [mm]	10
Sample port window material	BaF ₂
Output port diameter [mm]	12.7
Light emitting element	SiC globar
Source temperature [K]	~1'550
Rated life-time [hours]	10'000
Power consumption (7-12 VDC) [W]	30
External dimensions [mm]	93x93x130
Weight [g]	1300

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
 Please contact info@arcoptix.com for more information.

BATTERY OPERATED PORTABLE FTNIR SETUP



Rugged backpack and handheld fibered probe




Embedded software for data acquisition

The ARCOptix FTNIR-PORT is the ideal platform for measuring diffuse reflectance (DR) in the field. Battery operated, it features a light source, spectrometer, fibered reflection probe as well as an active cooling system with dust filters for robustness and long lasting operation. The entire data acquisition and processing is performed on the fly, and can be conveniently operated from a dedicated handheld measuring head. Easy data traceability is ensured thanks to the addition of metadata (GPS coordinates). The fibered probe is mechanically resistant owing to its stainless steel protective sleeve and contact ferrule. The casing can be easily opened to access the battery in order to swap it or recharge it.

Specifications

Product code	FTNIR-PORT
Instrument	FTNIR-L1-025-2TE (OEM version)
Light source	ARCLIGHT-NIR (OEM VERSION)
Battery autonomy	3h (in operation)
Fibered probe	6x low-OH, 600 µm core diameter silica fibers, stainless steel protection sleeve
Measuring head	Handheld, acquisition trigger, contact ferrule
Backpack dimensions [mm]	450x320x170
Total weight [kg]	4.5
Software interface	Embedded with GPS coordinates and personalized file saving system
Cooling system	Active (air flow) with dust filters

FIBERED REFLECTION PROBES




Fibered reflection probes are common accessories to measure specular or diffuse reflection in FT-IR systems. They are especially well suited for applications requiring a lot of flexibility. We offer fibered probes of three different kinds, depending on the fiber constituting glass, each targeting a specific spectral range. Light is brought to the sample via a peripheral fiber bundle which helps in achieving homogeneous illumination, and is collected by single, central fiber.

Specifications

Product code	R7-NIR-600-200F	R2-IFG-500-150F	R7-PIR-900-150F
Spectral range [μm]	0.4-2.5	2-5.5	4-16
Spectra range [cm⁻¹]	4'000-25'000	1'800-5'000	650-2'500
Illuminating fiber bundle	6x low-OH, 600 μm core diameter silica fibers	2x 500 μm core diameter IFG fibers	7x 400 μm core diameter PIR fibers
Reading fiber	low-OH, 600 μm core diameter silica fiber	500 μm core diameter IFG fiber	900 μm core diameter PIR fiber
Probe ferrule	Stainless steel ferrule ø1/8" (3.175mm x 74mm)	Stainless steel ferrule ø1/4" (6.35mm x 74mm)	Stainless steel ferrule ø1/4" (6.35mm x 74mm)
Length [m]	2	1.5	
Fiber connectors	SMA905		
Recommended Instrument	FTNIR-L1-025-2TE	FTMIR-L1-060-4TE	FTMIR-L1-160-LN2

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
 Please contact info@arcoptix.com for more information.

OPTICAL FIBERS



Optical fibers are a common building block of many optical systems. Low-OH fused silica multimode fibers are available for operation in the near-infrared. In the mid-infrared, we offer either IFG fibers or polycrystalline fibers depending on the spectral range of interest. Our fibers are terminated with SMA905 connectors and are protected by a polymer jacket. Other connector types, protection jacket or lengths are available on request.

Specifications

Product code	FIB-NIR-600-100	FIB-IFG-600-100	FIB-PIR-900-100
Spectral range [μm]	0.4-2.5	2.0-5.5	4-16
Spectra range [cm ⁻¹]	4'000-25'000	1'600-5'000	650-2'500
Fiber material	Low-OH fused silica	Indium Fluoride	Polycrystalline
Core diameter [μm]	600	600	900
Glass refractive index	1.44	1.5	2.15
NA	0.22	0.3	0.25
Maximum operating temperature [°C]	125	100	140
Minimum bending radius [mm]	100	100	130
Length [m]	1*	1*	1*

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE

* other lengths available on request

Please contact info@arcoptix.com for more information.


ARCOPTIX CUV



The ARCOptix CUV is our homemade cuvette holder system for transmission measurements in the near infrared. It features two fiber collimators with SMA905 connectors, a 12.5 mm wide sample port that accepts cuvette with an optical path up to 10 mm long, as well as a filter slot to avoid spurious effects (heating of the sample, extra shot noise) resulting from a too high illumination.

Specifications

Product code	PART-CUV-NIR
Spectral range [μm]	0.4-2.5
Spectra range [cm⁻¹]	4'000-25'000
Lens material	Borosilicate (BK7)
Fiber connector	SMA905





SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
 Please contact info@arcoptix.com for more information.

ARCOPTIX DRS98



The ARCOptix DRS98 is a reference element to be used in diffuse reflection measurements. Our DRS98 is made from optical grade Teflon (PTFE), which is most suited for the VIS-NIR range. Alternatively, our DRS98G, which is made from diffused gold, can be operated down to the mid-infrared.

Specifications

Product code	PART-DRS-NIR	PART-DRS-MIR
Reflectivity range (>90%) [μm]	0.4-2.5	0.5-16
Reflectivity range (>90%) [cm^{-1}]	4'000-25'000	625-20'000
Clear aperture [mm]	20	
Dimension [mm]	$\varnothing 38 \times 11$	
Reflecting material	 Optical grade PTFE	 Diffused Au

SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE
 Please contact info@arcoptix.com for more information.